

# ImmunoRun

Biogal Galed Labs.

## Canine /Feline Lateral Flow Test Kits For:

- ▶ FELINE IMMUNODEFICIENCY VIRUS Ab. (FIV)
- ▶ FELINE LEUKEMIA VIRUS Ag (FeLV)
- ▶ CANINE PARVO VIRUS Ag (CPV)
- ▶ CANINE DISTEMPER VIRUS Ag (CDV)
- ▶ CANINE HEARTWORM Ag (CDI)
- ▶ CANINE/FELINE GIARDIA Ag (CFG)
- ▶ CANINE LEISHMANIA Ab (CLI)

### Rapid

- One step procedure with accurate results within 10 Minutes

### Reliable

- High sensitivity, high specificity, high accuracy

### Easy


- Easy to use and easy to read results

### and Flexible

- Economical and inexpensive test-kit, suitable for the average sized clinic or larger clinic/lab - convenient pack sizes of 5 or 50 tests.



- \* No Refrigeration needed – Easy to send and store, no waiting for kit to warm
- \* Long shelf life of 24 months



Biogal Galed Labs.

# ImmunoRun


Biogal Galed Labs.

## ImmunoRun Performance Data Sheet - February 1, 2011

Name of Product	Sensitivity	Specificity	Comparison Test
CPV Ag	100%	98.8%	Hemagglutination Test (HT)
CHW Ag	*94.4%	100%	Necropsy
CDV Ag	100%	98.5%	RT-PCR
FIV Ab	96.8%	99.7%	Western blot
FeLV Ag	94.7%	99.7%	Virus isolation
Giardia Ag	100%	100%	Commercial Test
Leishmania Ab	95.6%	98%	IFA

\*CHW antigen sensitivity:

1 Female worm: 85%, 2 female worms: 100%, over 3 female worms:100%



**Biogal – Galed Labs.** Tel: 972-4-9898605 | Fax: 972-4-9898690

E-mail: [info@biogal.co.il](mailto:info@biogal.co.il) | Site: [www.biogal.co.il](http://www.biogal.co.il)