



R-CU510

Versatile Animal ICU

RWD



Quiet, Comfortable & Safe



The reference data can be calculated by selecting animal type-age-disease name



Touchscreen LCD display for best simplicity



With heating and cooling capabilities



Automatic ventilation based on CO₂ concentration parameters



Precise control of the O₂ concentration



Air filter • UVC lamp • Anion function for effective sanitation management, to avoid patients cross infection



Tempered glass windows for better observation



Built-in failure alarm covering abnormal temperature, humidity and CO₂ levels.

RWD Vet Anesthesia, Easy & Safe

R-CU510 Versatile Animal ICU

Versatile Animal ICU is indicated for postoperative recovery, shocks, head injury, neurologic disease, tissue hypoxia, cardiopulmonary disease, and numerous other symptoms and situations.

It provides animal patients with the most comfortable medical care environment to help animals recover faster.



Specifications

Item	Parameters
Temperature control range	15-40°C(68~104°F), accuracy :±0.5°C
Humidity control range	30~70%RH, accuracy :±1%RH
Ventilation	Level 1~Level 6
LED lamp	Level 1~Level 3
UVC lamp	5~60min(adjustable)
CO2 concentration monitoring	Range :500-4000ppm measurement accuracy : ±10ppm
O2 concentration control range	20~60%
IR lamp	Power:100W, 1~5min(adjustable)
Anion	7.2 million/cm3, 5~30min(adjustable)
Noise	≤50dB
Size	External: W930*H790*L536mm Internal: W690*H490*L470mm
Power supply voltage	AC 100~120V / 220~240V, 50/60Hz

RWD Life Science Co.,Ltd

Add: 19-20/F, Building 7A, Shenzhen International Innovation Valley, Dashi 1st Road, Nanshan District, Shenzhen, Guangdong, P.R. China **Enquiries:** market@rwdstco.com

RWD Life Science Inc.

Add: 850 New Burton Road, Suite 201, Dover, DE 19904 (New Office)
Tel: (858)900-5879 **E-mail:** service@rwdls.com