

SUPERLAB[®]
Your lab – Our passion

mindray
animal medical

vetXpert C5
Veterinarski biohemijski analizator

Sveobuhvatna dijagnostika
Idealan izbor



Međunarodni standardi za preciznu detekciju

Smanjenje rizika od pogrešne dijagnoze



Međunarodno priznato za preciznu sledljivost

- Vrhunski proces proizvodnje reagensa za izuzetne performanse
- Uporedivo sa međunarodnim referentnim laboratorijama za precizne rezultate testiranja



Visoke performanse uz minimalne greške

- Više od 20 godina iskustva u razvoju biohemijskih instrumenata sa garantovanim kvalitetom
- Dinamička tehnologija standardizacije osigurava konzistentne performanse reagensa



Korisnički prilagođen dizajn omogućava brzo testiranje


- Jednostavno rukovanje bez kidanja, pritiskanja ili predzagrevanja
- Potrebno je samo 100 μ L uzorka za više testova, sa rezultatima za 12 minuta



Inteligentna automatska kontrola kvaliteta

- Opremljen sistemom za kontrolu kvaliteta u realnom vremenu 360° za bezbrižno korišćenje
- Praćenje kritičnih tačaka obezbeđuje trenutna upozorenja, pružajući vam potpuni mir

Generisanje 29 rezultata istovremeno uz minimalnu količinu uzorka



Test profili

Paneli	Parametri
Comprehensive Diagnosis Panel	TP, ALB, GLOB, A/G, ALT, AST, AST/ALT, ALP, GGT, TBIL, BUN, CREA, BUN/CREA, AMY, TC, GLU, Na ⁺ , K ⁺ , Na/K, Cl ⁻ , PHOS, Ca, CaxP, CK, LDH, TBA, TG, Mg, LIPA-D*
Diagnosis/Health Check Panel	TP, ALB, GLOB, A/G, ALT, ALP, GGT, TBIL, BUN, CREA, BUN/CREA, AMY, TC, GLU, PHOS, Ca, CaxP, CK, TG
Preanesthetic Evaluation Panel	TP, ALB, GLOB, A/G, ALT, ALP, BUN, CREA, BUN/CREA, GLU
Liver Recheck Panel	TP, ALB, GLOB, A/G, ALT, AST, AST/ALT, ALP, GGT, TBIL, BUN, TC, GLU, TBA
Kidney Recheck Panel	TP, ALB, GLOB, A/G, BUN, CREA, BUN/CREA, PHOS, Ca, CaxP, Na ⁺ , K ⁺ , Na/K, Cl ⁻
Electrolyte Panel	BUN, CREA, BUN/CREA, GLU, Na ⁺ , K ⁺ , Na/K, PHOS, Ca, CaxP, Cl ⁻ , tCO ₂ , Mg
Avian/Reptilian Panel	TP, ALB, GLOB, A/G, AST, BUN, TC, GLU, Na ⁺ , K ⁺ , Na/K, PHOS, Ca, CaxP, CK, LDH, UA, Mg, TBA, TG
Large Animal Panel*	TP, ALB, GLOB, A/G, AST, ALP, GGT, BUN, PHOS, Ca, CaxP, CK, Mg
Equine Panel*	TP, ALB, GLOB, A/G, AST, ALP, GGT, TBIL, BUN, CREA, BUN/CREA, GLU, Na ⁺ , K ⁺ , Na/K, Ca, tCO ₂ , CK, LDH
Canine Inflammation Panel*	BUN, CREA, BUN/CREA, AMY, c-CRP*, LIPA-D*
Feline Inflammation Panel*	ALB, ALP, GGT, BUN, CREA, BUN/CREA, TBA, f-SAA*, LIPA-D*

* Uskoro

Specifikacija

Tip uzorka	plazma, serum, puna krv
Zapremina uzorka	100uL
Vreme do rezultata	≤12min
Ekran	7-inčni podesivi ekran osetljiv na dodir
Dimenzije	195 × 295 × 250 (mm)
Težina	≤ 5kg

Nadogradivo imuno modulom u all-in-one uređaju.

Podrška za povezivanje sa vetXpert sistemom za upravljanje.



LinkedIn Facebook Instagram

Search Mindray Animal Medical


www.mindrayanimal.com
service@mindrayanimal.com

P/N:ENG-vetXpert C5-210285X4P-20250214

©2025 Shenzhen Mindray Animal Medical Technology Co., LTD. All rights reserved.

SUPERLAB[®]
Your lab - Our passion

Milutina Milankovića 25,
11070 Novi Beograd
Tel/Fax: +381 11 22 22 222
Mob: +381 64 8720 137 / 146, 199, 210
E-mail: vetlab@super-lab.com
www.super-lab.com
www.vetlabprodaja.com



mindray
animal medical